

# Large Camera for DESY FLASH XFEL Studies

#### Introduction

Xcam specialises in producing custom and prototype CCD cameras for leading-edge science experiments all over the world. In June 2008, XCAM was asked to design and produce a large area CCD camera for FLASH beam-time in October 2008.



3 Camera units shown in the configuration that was designed for the DESY FLASH experiments

Figure 1:

Total number of 15 micron pixels for the dual camera is 4096 x 4096

Total image area: 61.4mm x 61.4mm

## **Multiple Camera Combinations**

The design of the camera permits use in a number of configurations, such as the one shown above which is suitable for SAXS/WAXS type diffraction studies, where both small and wide angle diffraction detail needs to be captured.



Figure 2

(Above) Single unit Camera and (right) XCAM single unit on base plate in foreground with sensor exposed to view



## **Extended Dynamic Range**

Systems offered provide novel techniques that offer dynamic range extension resulting in *single photon detection* in outer part of image combined with *full-well capability* in centre regions of image.

#### **Slit Separation and Adjustment**

The front dual camera is on a base plate with an adjustable slit: minimum slit width is 300 microns, maximum slit width is user-definable Distance from edge of image area to edge of CCD package is 300 microns for each CCD



Figure 3: Showing dual camera unit being commissioned with (above) blue finger plates which protect CCDs until required, and (below) single unit camera being tested individually



Direct Detection from 100eV to 20keV

**Multi-Synchronisation sequencer cards** permit multiple CCDs to be operated in synchronisation for low noise performance.

Multiple CCD camera systems operated with a single PC through XCAM software., or capability for users to develop their own software calling XCAM dlls.

**Software/Hardware Triggering** to start Erase and Integration sequences, for synchronisation with experiment. Custom triggering schemes .

See www. For presentation of results obtained by Henry Chapman's group with XCAM camera at FLASH see http://www.newlightsource.org/cambridge.html

XCAM Limited, Grove Farm, Moulton, Northampton, NN3 7TG, UK

Tel: 44 (0)1604 670729, Fax: 44 (0)1604 671570 Web: www.xcam.co.uk, Email: info@xcam.co.uk