

# Never compromise when your safety is at stake

Exclusively available in the UK from Helapet, the Berner range of cytotoxic resistant PPE and spill control products offer the highest level of cytotoxic protection against commonly used drugs within the healthcare industry. Reducing the risk of occupational exposure, the Berner PPE and Spill Kit range safeguards the user using type-tested barrier protection.

Berner protective clothing has been designed and stringently tested against European standards, providing an attractive combination of performance, breathability and comfort for professionals working closely with cytotoxic substances. Clothing options include gloves, gowns and oversleeves, each manufactured, packaged and sterilised for use in cleanroom environments.

Specialist Berner Cytotoxic Spill Kits and spill control products are capable of protecting the user against the risk of hazardous spills or leakages in controlled working conditions. Using highly absorbent materials in combination with a selection of cytotoxic resistant components, the Berner spill control range delivers optimal safety at a cost effective price.

All Berner products are supported by full certification and testing data with free samples available to order from Helapet on selected items.

#### Index

Latex Protective Gloves	3
Neoprene Protective Gloves	4
Nitrile Protective Gloves	5
Protective Gown	6
Protective Oversleeves	7
Cytotoxic Spill Kit	8
Berner Cytotoxic Spill Kit XP	9
Absorbent Workmat	10
ChemoSorb Pad	10

Chemical Permeation Table 11







#### **Features and Benefits**

- Purpose designed for handling cytotoxics
- Over 8 hours of permeation resistance\*
- Strong, durable design
- Fully textured for good grip of products
- Ready folded long 300mm cuff
- Sterile double bagged poly packaging facilitating smooth cleanroom transfer

\*permeation data page 9

# Latex Protective Gloves

#### Comfortable yet safe

Berner latex gloves stand out due to their exceptional elasticity, high wearable comfort and top protective properties against cytotoxics, chemicals and biological agents. Each glove is textured for superior grip, with a long beaded cuff for added protection above the wrist.

#### Material

- Natural latex (powder-free)
- Colour: blue
- Very low protein: 17 μg/gHypoallergenic: 0.5 μg/g
- Polymer coating for easy donning and double gloving

#### Compliance

- Category III PPE
- Type-tested and certified
- Comprehensive permeation tests (as per EN 374-3)
- Viral resistance test: ASTM F1671
- AQL= 1.0 (Water test)
- 5 year shelf life
- Ideal for GMP compliant production of cytotoxics

#### Ordering Information

GL3010	GL3025
6.5 (S) • sterile • 20 x 10 pairs/pck	8 (ML) • sterile • 20 x 10 pairs/pck
GL3015	GL3030
7 (SM) • sterile • 20 x 10 pairs/pck	8.5 (L) • sterile • 20 x 10 pairs/pck
GL3020	GL3040
7.5 (M) • sterile • 20 x 10 pairs/pck	9 (XL) • sterile • 20 x 10 pairs/pck

Try before you buy and request **free samples** of our Berner cleanroom gloves. Order your preferred glove and size online or email sales@helapet.co.uk



#### **Features and Benefits**

- Purpose designed for handling cytotoxics
- Over 8 hours of chemical resistance\*
- 100% latex free
- Tear resistant
- Ready folded for easy donning and prevent contamination
- Sterile double bagged poly packaging facilitating smooth cleanroom transfer

## Neoprene Protective Gloves

## Where latex sensitivity is an issue

Neoprene protective gloves are made from high quality polychloroprene and offer excellent permeation resistance against cytotoxics, chemicals, biological agents and viruses. Owing to the material composition, these gloves are particularly suitable for people suffering from latex allergies, without compromising on comfort or protection performance.

#### Material

- Polychloroprene
- Colour: Beige
- Powder-free
- No detectable allergens
- Polymer coating for easy donning

#### Compliance

- Category III PPE
- Type-tested and certified
- Comprehensive permeations tests (as per EN 374-3)
- Viral resistance test: ASTM F1671
- AQL= 1.0 (Water test)
- 5 year shelf life
- Ideal for GMP compliant production of cytotoxics

#### Ordering Information

GP3060100	GP3080100
6.5 (S) • sterile • 10 x 10 pairs/pck	8 (ML) • sterile • 10 x 10 pairs/pck
GP3070100	GP3085100
7 (SM) • sterile • 10 x 10 pairs/pck	8.5 (L) • sterile • 10 x 10 pairs/pck
GP3075100	GP3090100
7.5 (M) • sterile • 10 x 10 pairs/pck	9 (XL) • sterile • 10 x 10 pairs/pck



<sup>\*</sup>permeation data page 9

# Nitrile Protective Gloves

## Cytotoxic resistance and sensitivity combined

Setting the standard in comfort and exceptional glove protection, the Berner nitrile gloves are the ideal choice for long work periods with hazardous substances.

Manufactured from a copolymer of latex free materials grants over 8 hours of validated cytotoxic resistance, combined with outstanding tactility.

#### Material

- Special nitrile mixture
- · Colour: Light Blue
- No detectable allergens
- · Cuff length: 300mm

#### Compliance

- Tested with 15 cytotoxic substances
- Category III PPE
- Comprehensive permeation tests (as per DIN EN 374-3)
- Viral resistance test: ASTM F1671
- AQL= 1.5 (Water test)
- 3 year shelf life
- Terminally gamma sterilised

#### Ordering Information

#### Non sterile

#### GN3010

6 (XS) • non sterile • 100 pcs/pck

GN3015

7 (S) • non sterile • 100 pcs/pck

GN3020

8 (M) • non sterile • 100 pcs/pck

GN3025

9 (L) • non sterile • 100 pcs/pck

GN3030

10 (XL) • non sterile • 100 pcs/pck

#### Sterile

#### GN9610100

6 (XS) • sterile • 10 x 10 pairs/pck

GN9615100

7 (S) • sterile • 10 x 10 pairs/pck

GN9620100

8 (M) • sterile • 10 x 10 pairs/pck

GN9625100

9 (L) • sterile • 10 x 10 pairs/pck

GN9630100

10 (XL) • sterile • 10 x 10 pairs/pck



#### **Features and Benefits**

- Over 8 hours of chemical resistance\*
- Optimal handling with textured finger grips
- Ambidextrous, latex free presentation
- Long 300mm beaded cuff for secure fit
- Sterile and non sterile for different work applications

\*permeation data page 9





### **Protective Gown**

#### Proven front and arm protection

Lightweight, low linting and polycoated, the Berner Protective Gown grants frontal bodily protection against cytotoxics and biological agents in controlled working environments. Sealed seams deliver optimum protection, with knitted cuffs and a breathable back for additional comfort. Individually packed and sterilised for use in controlled environments, each gown can be worn over cleanroom coveralls.

#### Material

- Polypropylene spun bond (ppsb)
- Liquid impermeable coating on front and sleeves
- Colour: Blue
- I atex free

#### Compliance

- Permeation tested in accordance with EN 374-3:2003
- Category III PPE
- 4 year shelf life
- Terminally sterilised

#### Ordering Information

AB2010

AB2030

Small • sterile • 10 pcs/pck

Large • sterile • 10 pcs/pck

AB2020

Medium • sterile • 10 pcs/pck

AB2040 X Large • sterile • 10 pcs/pck

#### **Features and Benefits**

- Assured protection against cytotoxic agents\*
- Lint-free material
- Secure Velcro® fastening at the neck
- Double sewn wrapping bands
- · Comfortable knitted cuffs
- Sterile for cleanroom use

#### Size Chart (cm)

Size	Chest	Neck	Sleeve Length	Body Length	Tie Distance from Hem	Size Around
S	71.1	53.3	80	104	68.6	127
М	71.1	58.4	80	114	68.6	127
L	76.2	66	83.8	114	68.6	152
XL	83.8	68.6	88.9	117	68.6	167.6



<sup>\*</sup>permeation data page 9

## Protective Oversleeves

# Specific arm protection from cytotoxics

Granting maximum arm protection against low and high risk cytotoxic substances, Berner Oversleeves offer excellent permeation breakthrough times against a list of cytotoxic and chemical agents. Terminally sterilised and shaped for universal fit over cleanroom coveralls, makes them ideal for preparation activities in controlled working environments.

#### Material

- Polypropylene spun bond (ppsb) with poly-coating
- Ultrasonically sealed seams
- Colour: BlueLength: 52cmLatex free

#### Compliance

- Permeation tested in accordance with EN 374-3:2003.
- Category III PPE
- 4 year shelf life
- Terminally sterilised

#### Ordering Information

#### AB0001

52cm • sterile • 40 pairs/pck



#### **Features and Benefits**

- Impermeable arm protection against cytotoxic agents\*
- Liquid proof coating
- Lint-free material
- Comfortable knitted cuffs
- Universal fit
- Sterile for cleanroom use

\*permeation data page 9



## Cytotoxic Spill Kit

# The rapid, safe solution for dealing with cytotoxic spills

The ideal grab and go emergency spill kit containing a quality selection of compliant Personal Protective Equipment (PPE) and necessary aids to contain and dispose of hazardous cytotoxic spills. All contents are packed in the order the user will need them, enabling quick and safe access. Each spill kit comes with

instructions for use and an illustrative guide. Helapet offer four convenient spill kit configurations to address glove and gown size preferences.

#### Ordering Information

#### CS4004

Includes:

size 7.5 latex glove and medium gown non sterile • 2 pcs/pck

#### CS4004N

Includes:

size 7.5 neoprene glove and large gown non sterile • 2 pcs/pck

#### CS4004A

Includes:

size 9 latex glove and large gown non sterile • 2 pcs/pck

#### CS4004L

Includes:

size 9 neoprene glove and large gown non sterile • 2 pcs/pck

#### **Features and Benefits**

zic Drug SPIL

Zytostatika SPILL KIT

Cytotoxic Drug SPILL KIT

- Contains list of essential protective garments and accessories
- Validated permeation resistant PPE
- Enclosed accident report and illustrative guide
- Four gown and glove size combinations to suit

#### Contains

- Cytotoxic protective gown
- Cytotoxic protective gloves
- Protective gloves
- Protective overshoes
- FFP3-respiratory protection mask
- Protective goggles

- 1 x ChemoSorb pad
- 6 x Isysoft cloths
- Waste removal bag
- Hazardous waste removal bag, white
- 2 x cable ties
- User information

### Free Training Sessions

Helapet provide free on-site training sessions, introducing you and your staff to the benefits and proper use of our Spill Kit products, enquire today to arrange your dedicated training session.



Spilt cytostatic drugs!

# Cytotoxic Spill Kit XP

Maximum safety protection for the disposal of wet and dry cytotoxic spills

The Berner XP Spill Kit perfectly equips you with everything you need to safely contain and rapidly clean up accidental cytotoxic spills, using completely latex free components to reduce the risk of allergic reactions.

All contents are packed in the order you

require them for quick and easy accessibility, with supporting documentation on the correct use of the spill kit.



#### CS4124

medium • non sterile • 1 pc/pck

#### CS4133

large • non sterile • 1 pc/pck

#### CS4135

XX large • non sterile • 1 pc/pck

#### Contains

- Cytotoxic protective overall
- Cytotoxic protective gloves
- Protective overgloves
- Cytotoxic protective overboots
- FFP3 respirator mask
- Protective goggles
- Marker
- Warning sign
- 1 x ChemoSorb pad

- 6 x Isysoft cloths
- Scoop
- Sweeper
- Tongs
- Distilled water
- Waste bag
- Special waste bag, white
- Cable ties
- User information

#### **Features and Benefits**

• Suitable for wet and dry spills

#### Free Spill Kit Poster

Request your free 11 point procedure poster on the proper emergency application of Berner Cytotoxic Spill Kits.







#### **Features and Benefits**

- Lint free 3 layer protective system
- Highly absorbent
- Non-slip backing
- Reliable cytotoxic barrier properties
- Latex free

# Absorbent Workmat

# Safe work surface protection for cytotoxic drug preparations

Specially designed for use with cytotoxic preparations, the 3 ply blue mesh structure provides effective absorption and containment of accidental drips and spills. Lint free, sterile packaged and available in two convenient sizes, the Berner Absorbent Workmat is ideal for use in cleanroom and isolator environments.

#### Ordering Information

#### WP8001

26 x 42cm sterile • 50 pcs/pck

#### WP8101

41 x 56cm sterile • 40 pcs/pck



#### **Features and Benefits**

- Absorbs up to 3L of liquid
- Suitable for cytotoxic drugs
- Low linting
- Latex free
- Easy to handle

### ChemoSorb Pad

## Extra absorption when and where you need it

An economical solution for the safe absorption and containment of liquid spills. Each ChemoSorb Pad can reliably absorb chemo and non-chemo related spills, binding with the liquid to create a gel like mass which allows for easy disposal.

#### Ordering Information

#### CS120016

38cm x 54cm non sterile • 10 pcs/pck



### Permeation Breakthrough Times (min)

Chemical	Berner Latex Glove	Berner Neoprene Glove	Berner Nitrile Glove	Berner Cytotoxic Gown	Berner Oversleeves
5-Fluorouracil	>60	>480	>480	> 480	> 480
Acetone		20			
Acetonenitrile		20			
Amsacrin				> 480	> 480
Benzine	>180		45		
Bleomycins					
Butanone (MEK)		2			
Carboplatin	>180	_			
Carmustine	>60	90	>480	> 140	> 140
Chlorhexidine (CHX)	>00	30	>480	> 140	> 140
Cisplatin	>60	>480	>480	> 480	> 480
Cycloposphamide	>240	>400	>480	> 480	> 480
Cycloposphamide monohydrate	>60	>480	>400	> 400	<i>&gt;</i> 400
	>00	>400	× 400		
Cytarabine	. 040		>480		
Dacarbazine	>240		. 400		
Daunorubicin hydrochloride	>60		>480		
Denatured Alcohol	00	45			
Diethylamine	>30	45			
Doxorubicin hydrochloride	>180	>480		> 480	> 480
Ethanol		2	26		
Ethidium bormide (1%)			>480		
Etoposide	>240		>480		
Formaldehyde (10%)					
Formaldehyde (30%)					
Formaldehyde (35%)		>480	>480		
Glutaraldehyde (1%)	>480	>480	>480		
Glutaraldehyde (2.5%)					
Haxane		<7	40		
Hydrogen peroxide (30%)			>480		
Hypochloric acid (32%)					
Ifosfamide					
Irinotecan hydrochloride			>480		
Isopropanol	>30	>35	48		
Isopropanol/Carmustine	>30 / >120				
Methanol	7 00 / 7 120	20			
Methotrexate	>60	20	>480	> 480	> 480
Methotrexate + Ethanol	>00		> 100	> 400	> 100
Methyl medthacrylate (MMA)		<2			
	>180	ζζ	>480		
Mitomycin	>100		>400	<b>\ 100</b>	< 100
NaOH 30%			~ 100	> 480	> 480
Oxaliplatine			>480		
Povidone Iodine (10%)	0.40		400	400	400
Paclitaxel	>240		>480	> 480	> 480
Perchloic Acid (70%)		100	>480		
Sodium hydroxide (10%)	>480	>480			
Sodium hydroxide (30%)	>60				
Sodium hydroxide (40%)	>60	75			
Sodium Hypochlorite (12%)					
Sulphuric Acid (40%)	>480	>480			
Sulphuric Acid (96%)	>30	45	14		
Thiotepa	29		>480	> 480	> 480
Vinblastine	>180				
Vincristine	>180	>180	>480	> 480	> 480
Xylene		<2			



#### **Helapet Limited.**

Lister House, Blackburn Road, Houghton Regis, LU5 5BQ, UK T: +44 (0)1582 501980 F: +44 (0)1582 501981 E: sales@helapet.co.uk

helapet.co.uk/berner

#### Follow us







