

RE-USABLE ELECTROSURGERY

CATALOGUE

















We design and manufacture

Monopolar



Universal Fingerswitch

Our new range of fingerswitches which fit 2.4 and 4 mm electrodes

Page 2



Eco-switch

A unique concept in environmentally friendly semi-disposable fingerswitches which fit 2.4 and 4mm electrodes.

Pages 3-4



60+, 100+ & 200+ Fingerswitches

Page 5



Footswitch controlled handles

Footswitch controlled handles which fit 2.4 and 4mm electrodes

Page 6



Adaptors

Adaptors for all fingerswitches

Page 7



Universal Monopolar Cable

One cable that fits most instruments!

Pages 17-18



Monopolar Cables

A full range of monopolar cables to fit all instruments and generators

Pages 19-20



AUTOClabels

Autoclave indicator labels

Pages 27-28



Dispersive cables

Page 29



Monopolar forceps

Pages 31-32

Electrodes



Blades

Page 9



Needles

Page 10



Loops


Page 11



Balls

Page 12

electrosurgery accessories!



Bipolar





















Universal fingerswitch






Ergonomic, non-slip grips on the underside



Tactile, double insulated switches



Colour coded **AUTOClabel** tag is used to monitor usage with the permanently attached AUTOClabel punch.



Ordering Details:

Pack of 5 universal fingerswitches with blade electrodes, AUTOClabel and punch. Each fingerswitch is autoclavable sixty times:

PMS 2055 with 3m cable

PMS 2056 with 5m cable



Terminated with connector to fit:-Valleylab, BOWA Erbe ICC International Erbe VIO. Eschmann Conmed

PMS 2057 with 3m cable

PMS 2058 with 5m cable



Terminated with connector to fit:-4mm (female) Aesculap, Berchtold Erbe T Series, Huttinger, Martin, Micromed, Siemens Sutter BM780

PMS 2061 with 3m cable

PMS 2062 with 5m cable



Terminated with connector to fit:-4mm Erbe T-series

PMS 2059 with 3m cable

PMS 2060 with 5m cable



Terminated with connector to fit:-5mm ACC (standard) ICC, Erbe V10* Karl Storz *Depends on generator socket

Patented design











2.38mm and 4.00mm

(USA and European fittings)

active electrodes








Eco-switch Environmentally-friendly electrosurgery





Eco-switch Environmentally-friendly electrosurgery

















60+, 100+ & 200+ fingerswitches

3m

PMS 550

200 Plus fingerswitch with blade electrode pack 1

plus

AUTOClabel!

PMS 750

100 Plus fingerswitch with blade electrode pack 5

plus

AUTOClabel:

PMS 950

60 Plus fingerswitch with blade electrode pack 5

piu

AUTOClabel

5m

PMS 550/5

200 Plus fingerswitch with blade electrode pack 1

plus

AUTOClabel

PMS 750/5

100 Plus fingerswitch with blade electrode pack 5

plus

AUTOClabel

PMS 950/5

60 Plus fingerswitch with blade electrode pack 5

plu

AUTOClabel:

200+

Autoclavable for a maximum of 200 cycles

100+

Autoclavable for a maximum of 100 cycles

60+



Autoclavable for a maximum of 60 cycles

AUTOClabels




Autoclave tracking system



for adaptors see page 7



Footswitch-controlled handles





Ergonomically designed for precise control



Non slip grips - for safety



AUTOClabels autoclave tracking system



Ordering Details:

Pack of 3 footswitch contolled handles with blade electrode, AUTOClabel and punch, terminated with electrosurgery generator connector:

PMS 651 with 3m cable

PMS 651/5 with 5m cable



connector to fit: 4mm Aesculap Berchtold Erbe T. Series Eschmann TD300, TD311, DS302

Terminated with

Sutter BM780

with 3m cable PMS 652/5 with 5m cable

PMS 652



connector to fit: 8mm Bard, Birtcher Bovie, Bowa Conmed, Eschmann Medtrex Valleylab

Terminated with

PMS 654 with 3m cable

PMS 654/5 with 5m cable



Terminated with connector to fit: 5mm ACC (standard)

ICC Erbe V10* Karl Storz *Depends on generator socket

PMS 653 with 3m cable

PMS 653/5 with 5m cable



Terminated with connector to fit: 4mm (female) Aesculap, Berchtold Erbe T Series, Huttinger Martin, Micromed Siemens Sutter BM780

















3-pin plug to Berchtold / Martin generators

PMS 800 Pack of 1 Fingerswitch adaptor



3-pin plug to Erbe ACC/ICC generators with 5mm pin

PMS 850 Pack of 1 Fingerswitch adaptor



3-pin plug to Erbe T series generators with 4mm pin

PMS 860 Pack of 1 Fingerswitch adaptor





Connector from fingerswitch







Safety guard: Any break in the outer black insulation is highlighted by the yellow insulation layer below



Safety Shield: A slip-on sleeve designed to prevent potential arcing between the electrode and fingerswitch.



Indicator grooves: Enables sterilisation records to be maintained using the appropriate systems.



Unique number: Enables sterilisation record to be made of the product. (up to 25 uses)















Blade electrodes





Needle electrodes



















Loop electrodes





Ball electrodes























Coloured (standard) fibre light cables





Fibre light cables

Available cable length

1.8 m Add 18 to code Add 23 to code 2.3 m 3.0 m Add 30 to code

Available fibre-light bundle diameter

Add 35 to code 3.5 mm Add 49 to code


Available fibre light cables colours

for Blue G for Green R for Red for Yellow Y GY for Grey

Generate your own fibre optic cable codes

The codes listed left are the first part of your product code. To generate the second part, you need to know the length and diameter of the cable you require (see above), finally add the colour you wish (see above)

For example;



NOT FOR USE WITH XENON LIGHT SOURCES WITHOUT HEAT FILTERS


















Coloured high flex, armoured fibre light cables







Armoured cables

Crush resistant

- For use with Xenon light sources
- high-flex armoured cables
- choice of colours identify your department!

Fibre light cables

Available cable length

1.8 m	Add 18 to code
2.3 m	Add 23 to code
3.0 m	Add 30 to code

Available fibre-light bundle diameter

3.5 mm	Add 35 to code
4.9 mm	Add 49 to code


Available fibre light cables colours



Generate your own fibre optic cable codes

The codes listed left are the first part of your product code. To generate the second part, you need to know the length and diameter of the cable you require (see above), finally add the colour you wish (see above)


For example;



Product code: CPX2335/R

SUITABLE FOR USE WITH XENON LIGHT SOURCES



















Universal monopolar cables



Universal monopolar cables



Fits most male forcep ends



4.8mm Monopolar handles (foot controlled)



4.0mm Endoscopic



4.5 to 5.0mm Monopolar forceps



4.0mm Monopolar handle UK

AUTOClabels

Pop a hole in the AUTOClabels every time the product is autoclaved with the permanently attached AUTOClabels punch.

See page 28



















Monopolar cables

		Instrument connectors	Olympus resectoscope (4mm female instrument socket)	3mm endoscopic socket Olympus, Cook instruments	US screw fitting Olsen	
Genera	ator connectors		G	10	8	
4mm male	Aesculap Berchtold Erbe T. Series Eschmann TD300,		B-G 3m cable	B-10 3m cable	B-8 3m cable	
	TD311, DS302 Sutter BM780		B-G/5 5m cable	B-10/5 5m cable	B-8/5 5m cable	
8mm	Bard Birtcher Bovie Bowa		C-G 3m cable	C-10 3m cable	C-8 3m cable	
	Conmed Eschmann Medtrex and Valleylab		C-G/5 5m cable	C-10/5 5m cable	C-8/5 5m cable	
5mm	ACC (standard) ICC Erbe V10* Karl Storz		D-G 3m cable	D-10 3m cable	D-8 3m cable	
	*Depends on generator socket		D-G/5 5m cable	D-10/5 5m cable	D-8/5 5m cable	
4mm female	Aesculap Berchtold Erbe T Series Huttinger		G-G 3m cable	G-10 3m cable	G-8 3m cable	
	Martin Micromed Siemens and Sutter BM780		G-G/5 5m cable	G-10/5 5m cable	G-8/5 5m cable	



Monopolar cables

4mm socket endoscopic instrument	Adult Storz resectoscope	Comeg 2mm pin fitting Wolf Storz pin fitting	Childs Storz resectoscope	UK 4.8mm socket shrouded safety connector	Pop a hole in the A every time the procautoclaved with the attached AUTOCIA
7	6	4 (old 11)	3	1	
B-7 3m cable	B-6 3m cable	B-4 3m cable	B-3 3m cable	B-1 3m cable	
B-7/5 5m cable	B-6/5 5m cable	B-4/5 5m cable	B-3/5 5m cable	B-1/5 5m cable	
C-7 3m cable	C-6 3m cable	C-4 3m cable	C-3 3m cable	C-1 3m cable	
C-7/5 5m cable	C-6/5 5m cable	C-4/5 5m cable	C-3/5 5m cable	C-1/5 5m cable	
D-7 3m cable	D-6 3m cable	D-4 3m cable	D-3 3m cable	D-1 3m cable	
D-7/5 5m cable	D-6/5 5m cable	D-4/5 5m cable	D-3/5 5m cable	D-1/5 5m cable	
G-7 3m cable	G-6 3m cable	G-4 3m cable	G-3 3m cable	G-1 3m cable	These cables are for use with UK and
G-7/5 5m cable	G-6/5 5m cable	G-4/5 5m cable	G-3/5 5m cable	G-1/5 5m cable	Prima monopolar forceps Pages 31-32

labels

AUTOClabels oduct is he permanently Clabels punch.


















Universal bipolar cables





Universal bipolar cables

AUTOClabels

Pop a hole in the AUTOClabels every time the product is autoclaved with the permanently attached AUTOClabels punch.

See page 28



2mm Martin, Aesculap, Integra, Wolf, Berchtold, Sutter BM780	Older Eschmann 300 series (Lemo)	Twin 4mm pin, Erbe Int. Eschmann, Alsa MCD series	Aesculap (Older models GK50)	Storz	Bowa, EMC, Erbe Int. Ethicon, Lamidey, Valleylab Coa comp, Olympus, Soering	Codman CMC11	Valleylab (autobipolar models only)
	H	* *					
Q	R	S	Т	U	V (old N)	W	Υ
Q-30 4m cable	R-30 4m cable	S-30 4m cable	T-30 4m cable	U-30 4m cable	V-30 4m cable	W-30 4m cable	Y-30 4m cable















AUTOCIabels Autoclave tracking system PAGE 28	Generator connectors	Conmed	Dolley EMC	Erbe, Siemens Sutter BM50	2mm Martin, Aesculap, Wolf Berchtold, Sutter BM780
Instrument connectors		L	M	Р	Q
European fitting	12	L-12 3m cable	M-12 3m cable	P-12 3m cable	Q-12 3m cable
		L-12/5 5m cable	M-12/5 5m cable	P-12/5 5m cable	Q-12/5 5m cable
Karl Storz	14	L-14 3m cable	M-14 3m cable	P-14 3m cable	Q-14 3m cable
	14	L-14/5 5m cable	M-14/5 5m cable	P-14/5 5m cable	Q-14/5 5m cable
R Wolf	15	L-15 3m cable	M-15 3m cable	P-15 3m cable	Q-15 3m cable
	13	L-15/5 5m cable	M-15/5 5m cable	P-15/5 5m cable	Q-15/5 5m cable
USA twin pin (straight)	10	L-16 3m cable	M-16 3m cable	P-16 3m cable	Q-16 3m cable
•	16	L-16/5 5m cable	M-16/5 5m cable	P-16/5 5m cable	Q-16/5 5m cable
USA twin pin (angled)		L-16A 3m cable	M-16A 3m cable	P-16A 3m cable	Q-16A 3m cable
	16A	L-16A/5 5m cable	M-16A/5 5m cable	P-16A/5 5m cable	Q-16A/5 5m cable



Older Eschmann 300 series (Lemo)	Twin 4mm pin, Erbe Int. Eschmann, Alsa MCD series	Aesculap (Older models GK50)	Storz	Bowa, EMC, Erbe Int. Ethicon, Lamidey, Valleylab Coa comp, Olympus, Soering V (old N)	Codman CMC11	Valleylab (autobipolar models only)	Conmed (handswitching bipolar forceps)
R-12 3m cable R-12/5 5m cable	S-12 3m cable S-12/5 5m cable	T-12 3m cable T-12/5 5m cable	U-12 3m cable U-12/5 5m cable	V-12 3m cable V-12/5 5m cable	W-12 3m cable W-12/5 5m cable	Y-12 3m cable	N/A
R-14 3m cable R-14/5 5m cable	S-14 3m cable S-14/5 5m cable	T-14 3m cable T-14/5 5m cable	U-14 3m cable U-14/5 5m cable	V-14 3m cable V-14/5 5m cable	W-14 3m cable W-14/5 5m cable	N/A	N/A
R-15 3m cable R-15/5 5m cable	S-15 3m cable S-15/5 5m cable	T-15 3m cable T-15/5 5m cable	U-15 3m cable U-15/5 5m cable	V-15 3m cable V-15/5 5m cable	W-15 3m cable W-15/5 5m cable	N/A	N/A
R-16 3m cable R-16/5 5m cable	S-16 3m cable S-16/5 5m cable	T-16 3m cable T-16/5 5m cable	U-16 3m cable U-16/5 5m cable	V-16 3m cable V-16/5 5m cable	W-16 3m cable W-16/5 5m cable	Y-16 3m cable	N/A
R-16A 3m cable R-16A/5 5m cable	S-16A 3m cable S-16A/5 5m cable	T-16A 3m cable T-16A/5 5m cable	U-16A 3m cable U-16A/5 5m cable	V-16A 3m cable V-16A/5 5m cable	W-16A 3m cable W-16A/5 5m cable	Y-16A 3m cable	N/A















AUTOClabels Autoclave tracking	Generator connectors	Conmed	Dolley EMC	Erbe, Siemens Sutter BM50	Martin, Aesculap, Wolf Berchtold, Sutter BM780
system PAGE 28		4	-		Î
Instrument connectors		L	M	Р	Q
Aesculap (old style))	17	L-17 3m cable	M-17 3m cable	P-17 3m cable	Q-17 3m cable
		L-17/5 5m cable	M-17/5 5m cable	P-17/5 5m cable	Q-17/5 5m cable
Conmed handswitching forceps	18	N/A	N/A	N/A	N/A
Storz/Spectrum (small Fischer plug)	20	L-20 3m cable	M-20 3m cable	P-20 3m cable	Q-20 3m cable
	20	L-20/5 5m cable	M-20/5 5m cable	P-20/5 5m cable	Q-20/5 5m cable
Ethicon bipolar scissors	22	L-22 3m cable	M-22 3m cable	P-22 3m cable	Q-22 3m cable
	22	L-22/5 5m cable	M-22/5 5m cable	P-22/5 5m cable	Q-22/5 5m cable
Universal (USA & European fittings)	30 See page 22	L-30 4m cable	M-30 4m cable	P-30 4m cable	Q-30 4m cable



Older Eschmann 300 series (Lemo)	Twin 4mm pin, Erbe Int., Eschmann, Alsa MCD series	Aesculap (Older models GK50)	Storz	Bowa, EMC, Erbe Int. Ethicon, Lamidey, Valleylab Coa comp, Olympus, Soering	Codman CMC11	Valleylab (autobipolar models only)	Conmed (handswitching bipolar forceps)
R	S	T	U	V (old N)	W	Υ	J
R-17 3m cable R-17/5 5m cable	S-17 3m cable S-17/5 5m cable	T-17 3m cable T-17/5 5m cable	U-17 3m cable U-17/5 5m cable	V-17 3m cable V-17/5 5m cable	W-17 3m cable W-17/5 5m cable	N/A	N/A
N/A	N/A	N/A	N/A	N/A	N/A	N/A	J-18 3m cable J-18/5 5m cable
R-20 3m cable R-20/5 5m cable	S-20 3m cable S-20/5 5m cable	T-20 3m cable T-20/5 5m cable	U-20 3m cable U-20/5 5m cable	V-20 3m cable V-20/5 5m cable	W-20 3m cable W-20/5 5m cable	N/A	N/A
R-22 3m cable R-22/5 5m cable	S-22 3m cable S-22/5 5m cable	T-22 3m cable T-22/5 5m cable	U-22 3m cable U-22/5 5m cable	V-22 3m cable V-22/5 5m cable	W-22 3m cable W-22/5 5m cable	N/A	N/A
R-30 4m cable	S-30 4m cable	T-30 4m cable	U-30 4m cable	V-30 4m cable	W-30 4m cable	Y-30 4m cable	N/A

















AUTOClabels - autoclave tracking system





AUTOClabels - autoclave tracking system

Simple!

It is easy to see and monitor product usage

- Pop a hole in the AUTOClabel every time the product is autoclaved with the permanently attached AUTOClabel punch or
- Track your product's autoclave history by using the unique label number printed at the top of each AUTOClabeltag or
- Use the blank **AUTOClabel** to apply your own barcode.

The AUTOClabel punch & label are easy to attach to the product, no tools needed, just clip on

No tie needed - Just push and lock the punch on the cable and the label just feeds back into itself - takes seconds!














ACL₃₀

Autoclave cycles 30

Pack 100 Punch and label sets ACL40

Autoclave cycles 40

Pack 100 Punch and label sets ACL60

Autoclave cycles 60

Pack 100 Punch and label sets ACL100

Autoclave cycles 100

Pack 100 Punch and label sets

ACL000 Blank for barcode

















Dispersive cables

















Monopolar forceps

PMS 14-7007

Pack of 1 Serrated dissecting forceps Adson 100mm + tip 2mm Terminated with UK 4.8mm pin

PMS 14-7000

Pack of 1 Serrated forceps Adson 125mm + tip 2mm Terminated with UK 4.8mm pin

PMS 14-7012

Pack of 1 Serrated dissecting forceps McIndoe 150mm + tip 2mm Terminated with UK 4.8mm pin

PMS 14-7002

Pack of 1 Serrated forceps McIndoe 150mm + tip 2mm Terminated with UK 4.8mm pin

PMS 14-7022

Pack of 1 Serrated half teeth dissecting forceps 150mm + tip 2mm Terminated with UK 4.8mm pin

PMS 14-7014

Pack of 1 Serrated dissecting forceps 180mm + tip 2mm Terminated with UK 4.8mm pin

PMS 14-7004

Pack of 1 Serrated forceps Waugh 180mm + tip 2mm Terminated with UK 4.8mm pin





Monopolar forceps




















Bipolar forceps





Bipolar forceps







Bipolar forceps





Bipolar forceps













NonStick bipolar forceps

- Significant reduction of tissue adhesion to tips
- Reduced tissue trauma
- Easy to clean during use


The **NonStick** range of bipolar forceps offers significant advantages over standard forceps when used in soft tissue surgery.

The adhesion of "coagulated" tissue to the tips is substantially reduced due to the heat dispersion properties of the silver metal tipped forceps; silver is the best heat conducting material available and is the only material which is bio-compatble.


How the **NonStick** forceps eliminate adhesion:-

Standard forceps during


"Hot Spots" occur coagulation




"Hot Spots" start to occur during coagulation




The "Hot Spots" cause localised sticking



Heat is conducted away rapidly to prevent adhesion



Tissue seared to the forceps tips 'breaks' off



The forceps can be removed without tissue adhesion



NonStick bipolar forceps



120-030NS 0.3mm tip 120-100NS 1mm tip 120-200NS 2mm tip Pack 1 Straight forceps 150mm

160-030NS 0.3mm tip 160-100NS 1mm tip 160-200NS 2mm tip Pack 1 Bayonet forceps 180mm

180-030NS 0.3mm tip 180-100NS 1mm tip 180-200NS 2mm tip Pack 1 Bayonet forceps 200mm

140-030NS 0.3mm tip 140-100NS 1mm tip 140-200NS 2mm tip Pack 1 Straight forceps 200mm

142-030NS 0.3mm tip 142-100NS 1mm tip 142-200NS 2mm tip

Straight forceps – angled tip 200mm



















Twin Footswitches

- Non-slip operation
- Robust
- Easy to clean
- Non corrosive finish



Twin Footswitches

		9
Connectors to generators		Twin footswitch
	Conmed / Aspen	FS250
	Bard / Davol	FS041
	Valleylab	FS206
	Erbe	FS720
	Eschmann 300 & 400 series	FS207
	Eschmann 800 series	FS307
	Unterminated footswitches without connectors	FS355F
All cables 5m long		



FS3250

FS3041

FS3206

FS3720

FS3207

FS3307

FS3355

- Ideal for use in sitting position
- Non-corrosive plastic finish



Built in cable restraint



Integral carrying handle



Ideal for use with all theatre footwear

















Single Footswitches

Single Footswitches





Single Footswitches

		Yellow (cut)	Blue (coag)	Black	White
Connectors to generator	'S	single footswitch	single footswitch	single footswitch	single footswitch
	Valleylab			FS210	FS211
	Erbe				FS501
	Birtcher			FS555	
	Eschmann 300 & 400 series	FS109	FS117	FS718	FS719
	Eschmann 800 series	FS209	FS217		
	Unterminated footswitches without connectors	FS355C	FS355D	FS355E	FS176

All cables 3m long (longer lengths available on request)















Pneumatically-operated Footswitch



Non-slip base



Fully secure cable

FS2000

Footswitch for use with pneumatically operated generators

3 meter cable





Prima Medical Limited

Terms and conditions

All orders for the sale of goods are accepted by Prima Medical Limited on these terms and conditions, which supersede any previously published terms in a catalogue or elsewhere

Offers and acceptance

No obligation shall arise as a result of a catalogue entry, quotation or any other offer, until the company accepts the order from a customer, company or other organisation seeking to make a purchase.

Orders

Where possible, please place orders by mail or fax. If obliged to order by phone, please ensure you have all the relevant details to hand in order to speed up this process. All orders will be processed only if made with full payment, except for such organisations and companies that have established account facilities with Prima Medical Limited. Orders will be processed for account holders only on provision of an Official Purchase Order, which has been duly signed.

Prices

Prima Medical Limited reserves the right to alter prices at any time without prior notice. Prices charged will be as quoted when the order is received and accepted. If any prices have risen when an order is received, you will be contacted and invited to consider the new price, prior to processing your order. If the price has dropped, we will automatically process the sale at the new price.

Satisfaction guarantee

It is our policy that if you are not totally satisfied with a product within 7 days of receipt, you may return the goods. We only ask that the goods must be unused, in an 'as new' condition. They must be in their original packaging and any included literature should be enclosed

We reserve the right to apply a restocking charge of up to 25%. The guarantee does not apply to the sterile products, 'specials' or cut lengths of cable. No goods may be returned without a returns authorisation number, issued by our customer services helpline (+44 (0) 1453 752 626). Claims for goods damaged in transit must be made within 3 days of delivery. In the case of damaged items, always retain them for our representative or courier. None of these conditions affect your statutory rights.

Back Orders

We reserve the right to ship part orders, or to hold orders until all items are available as required. We shall only levy one delivery charge. We reserve the right to cancel any part of an outstanding order should securing supply become untenable.

Cancelling an Order

Orders can only be cancelled by prior agreement of Prima Medical Ltd. Where possible, and shipment from vendor or ourselves has not been made, we shall agree. Should a customer default on an order outside these terms we reserve the right to deduct the relevant costs of shipment and administration incurred by Prima Medical Ltd.

Returns

No goods can be returned to the company without its express consent. In the event of faulty or damaged goods being received by the purchaser, please contact and advise us of the problem, whereupon arrangements to rectify the situation will be made. Claims for products shipped faulty or damaged must be made within 3 working days.

Warranty

Prima Medical Ltd will, at its discretion, repair. replace or refund a product that develops a fault/s within the given guarantee period that is a result of faulty workmanship or design. In the unfortunate event of a fault arising, you should return the product to your nearest distributor for examination by the distributors' staff, or telephone the Order Line for advice. The warranty is void if any modification or attempted repair has been made to a product, or the product as been subjected to accidental damage. This warranty excludes any liability for consequential loss arising from defective products.

Law

Any business transacted by Prima Medical Ltd is governed solely by English Law and none of the conditions stated affect your statutory rights.

We accept all major credit and debit cards













In the interest of security please telephone our sales department to pay by credit or debit card.

Prima Medical Ltd

Unit 21 Daniels Industrial Estate Bath Road Stroud Gloucestershire GL5 3TJ United Kingdom

+44 (0) 1453 752626

+44 (0) 1453 756750

sales@prima-medical.com

www.prima-medical.com

Company names mentioned within this brochure are the property of their respective owners and are not to be taken as an endorsement of affiliation with Prima Medical Ltd. Specifications are subject to change without notice

All rights reserved. No part of this publication may be reproduced or transmitted in any form or by any other means, whether electronic or mechanical including by photocopy, recording or information storage and retrieval system, without permission in writing from the publisher. Authorisation to photocopy items for internal or personal use, or the internal or personal use of specific clients is granted by Prima Medical Ltd.

All Products are Latex Free Copyright © 2009, Prima Medical Ltd.













view our online catalogue at www.prima-medical.com

