

Prima
Medical Limited

ELECTROSURGERY
COLPOSCOPY



+44 (0) 1453 752626

+44 (0) 1453 756750

sales@prima-medical.com

www.prima-medical.com



Designed and Manufactured by Prima Medical Ltd

All rights reserved. No part of this publication may be reproduced or transmitted in any form or by any other means, whether electronic or mechanical including by photocopy, recording or information storage and retrieval system, without permission in writing from the publisher. Authorisation to photocopy items for internal or personal use, or the internal or personal use of specific clients is granted by Prima Medical Ltd.

All Products are Latex Free

Copyright © 2006, Prima Medical Ltd.

CONTENTS

Biopsy Punches & Gynae Hook	2
Loop Electrodes (Re-Usable)	3
Fingerswitches & Ball Electrodes (Re-Usable)	4
Loop Electrodes (Single Use)	5
Fingerswitches & Ball Electrodes (Single Use)	6



Prima Medical Ltd

Unit 21, Daniels Industrial Estate, Bath Road, Stroud
Gloucestershire, GL5 3TJ, United Kingdom



+44 (0) 1453 752626







+44 (0) 1453 756750



sales@prima-medical.com



www.prima-medical.com




CALL +44 (0)1453 752626 TO PLACE AN ORDER; OR VISIT OUR COMPLETE ONLINE CATALOGUE AT

www.prima-medical.com




BIOPSY PUNCH
 1 Biopsy Punch Tishler-Morgan Re-usable pack 1 PMS1000

GYNAE HOOK
 2 Gynae Hook 250mm Single Use pack 50 PMS1790



LOOPS




RE-USABLE COLPOSCOPY ELECTRODES

Size may vary by 5%

1	Loop Electrode, 8 mm × 8 mm loop, 11cm,	pack 5	PMS1509
2	Loop Electrode, 10 mm × 5 mm loop, 11cm,	pack 5	PMS1510
3	Loop Electrode, 10 mm × 10 mm loop, 11cm,	pack 5	PMS1511
4	Loop Electrode, 12 mm × 8 mm loop, 11cm,	pack 5	PMS1512

5	Loop Electrode, 12 mm × 12 mm loop, 11cm,	pack 5	PMS1501
6	Loop Electrode, 15 mm × 8 mm loop, 11cm,	pack 5	PMS1513
7	Loop Electrode, 15 mm × 15 mm loop, 11cm,	pack 5	PMS1502
8	Loop Electrode, 20 mm × 10 mm loop, 11cm,	pack 5	PMS1514

9	Loop Electrode, 20 mm × 15 mm loop, 11cm,	pack 5	PMS1503
10	Loop Electrode, 25 mm × 10 mm loop, 11cm,	pack 5	PMS1515
11	Loop Electrode, 25 mm × 15 mm loop, 11cm,	pack 5	PMS1517
12	Loop Electrode, 25 mm × 25 mm loop, 11cm,	pack 5	PMS1504



PRIMA SAFETY FEATURES Safety guard




Any break in the outer black insulation is highlighted by the yellow insulation layer below

PRIMA SAFETY FEATURES Safety shield




A slip-on sleeve designed to prevent potential arcing between the electrode and fingerswitch.

PRIMA SAFETY FEATURES Indicator grooves



Enables sterilisation record to be made on the product. (up to 25 uses)

UNIQUE NUMBER



Enables sterilisation records to be maintained using the appropriate systems.


FINGERSWITCH



1

Electrode not included

BALL ELECTRODES



2


Prima 00-1058
PMS 1408

3

Prima 00-1058
PMS 1420

1 2 3 4 5 6 7 8 9 10 11 cm

LOOP ELECTRODES





RE-USABLE FINGERSWITCH

1 100 Plus Fingerswitch, 3m cable terminated with 3 pin plug,
(Autoclavable for a maximum of 100 cycles) pack 5 PMS750

RE-USABLE BALL ELECTRODES

2 3 mm Ball Electrode, 11 cm, pack 3 PMS1408
3 5 mm Ball Electrode, 11 cm, pack 3 PMS1420




LOOPS





SINGLE USE COLPOSCOPY ELECTRODES

Size may vary by 5%

1 Loop Electrode, 10mm x 5mm loop, Pink, pack 25	PMS1740	5 Loop Electrode, 15mm x 15mm loop, Blue, pack 25	PMS1743	9 Loop Electrode, 20mm x 20mm loop, Beige, pack 25	PMS1744
2 Loop Electrode, 10mm x 8mm loop, Aqua, pack 25	PMS1750	6 Loop Electrode, 20mm x 12mm loop, Grey, pack 25	PMS1751	10 Loop Electrode, 5mm x 5mm square loop, Light Green, pack 25	PMS1735
3 Loop Electrode, 10 mm x 10mm loop, Green, pack 25	PMS1742	7 Loop Electrode, 20mm x 13mm loop, Yellow, pack 25	PMS1752	11 Loop Electrode, 10mm x 4mm square loop, Green, pack 25	PMS1731
4 Loop Electrode, 15mm x 8mm loop, Orange, pack 25	PMS1746	8 Loop Electrode, 20mm x 15mm loop, Red, pack 25	PMS1745	12 Loop Electrode, 10mm x 8mm square loop, Brown, pack 25	PMS1732






FINGERSWITCH



BALL ELECTRODES




SINGLE USE FINGERSWITCH

1 Fingerswitch with 3 pin plug, 3m cable, boxed 25 PMS255

SINGLE USE BALL ELECTRODES

2 Ball Electrode, 3 mm, 12cm pack 25 PMS1761

3 Ball Electrode, 5 mm, 12cm pack 25 PMS1762



view our online catalogue at

www.prima-medical.com

Call today
to receive your
FREE
Prima catalogue

