

► SampleSafe Lock Box with Insulation

available for purchase for your clinical trial.



Overview

The SampleSafe Lock Box is available with air insulation or foam insulation. The SampleSafe Lock Box provide superior sample protection. Included with the lock box is a convenient carrying handle and a keyed lock.

The superior thermal characteristics of the SampleSafe excel in maintaining sample storage temperatures in a variety of different ambient conditions. When used in combination with the optional cooling packs, sample temperatures within the SampleSafe can be safely maintained for up to 10 hours.

Operating Features

- Unmatched air or foam insulation properties guarantee the highest sample integrity possible
- UV and moisture resistant moulded design provides years of quality service
- Separate compartment for cold packs prevent unwanted contact with patient samples
- Folding carry handle simplifies transport inside or outside
- Recessed lock protected from ice build up
- Optional mounting bracket for securing the SampleSafe
- False bottom isolates cold packs from sample
- Optional accessories available to include: Replacement keyed lock, Replacement key, Combination lock, False bottom, Under door bracket kit, Anchor bracket kit, 6 foot anchor cable, Rubber feet kit and cold pack.

Specification

SampleSafe Lock Box - Air Insulated

- Dimensions: 13.5 h x 12 w x 12.5 d inches
34.29 hx 30.48 w x 31.75 d cm
- Weight: 6.7lbs / 2.92 Kg

SampleSafe Lock Box - Foam Insulated

- Dimensions: 13.5 h x 12 w x 12.5 d inches
34.29 hx 30.48 w x 31.75 d cm
- Weight: 8.5lbs / 3.86 Kg



Please contact us for a quotation for the purchase of this product.

Product description - SampleSafe Lock Box



UK Office

T: +44(0)8456 777001

F: +44(0)8456 777002

E: sales@woodleyequipment.com www.woodleyequipment.com

USA Office

T: 1 800 471 9200 / 1 508 625 1693

F: 1 508 625 1721

WOODLEY
EQUIPMENT COMPANY LTD.